

## Reading data from WCF files

The purpose of this document is to outline the structure of the .WCF file format. The recommended method of accessing data from a WCF file is to use the DataRay OCX to open the file and retrieve the calculated values as outlined in a number of tutorials [on our website](#). However, some customers need to be able to access the WCF files byte-wise.

### File Structure

A WCF file contains 2 structures. First, there is a 5592 byte at the top of the file WC\_IMAGE\_DATA\_HEADER\_2. The first DWORD in this header contains the characters DRI, which should be used to check if the file is valid. After this header, there is a new structure for each frame stored in the .WCF file WC\_IMAGE\_DATA. This structure contains 944 bytes of header information before the actual image data. The image data is stored as a 1 dimensional array of 2-byte values. They are ordered by row. You will need to know the size of the data (rows x columns). These sizes can be read in from the Width and Height variables respectively.

```
typedef struct{
    DWORD Signature; //”DRI.”
    DWORD Type;
    DWORD Size;
    DWORD Images;
    DWORD ImagesSize;
    char Version[40];
    DRI_SETTINGS Settings;
} WC_IMAGE_DATA_HEADER_2;

typedef struct{
    int Signature;
    int Type;
    int Index;
    int Beams;
    int Size;
    int Width;// Number of horizontal pixels
    int Height;// Number of vertical pixels
    int CameraUpdateNumber;
    double XpixelSize; //Pixel horizontal
    double YpixelSize; //Pixel vertical
    int Bits;//Normal = 16
    int Key;
    int Peak;
    int Xoffset;// x start offset (unused pixels)
    int Yoffset;// y start offset (unused pixels)
    int Xlimit;// imagers total number of x pixels
    int Ylimit;// imagers total number of y pixels
    int OreintationDone;
    CPoint pPeakCenter;
    double DefinedFluencePower;
    double pUserCentroid[2];
    double Centroid[2];
    double GeoCentroid[2];
    double Baseline;
    double UserCentroid[2];
}
```

```

double      GeoCenter[2];
double      PeakCentroid[2];
double      Orientation;
double      Ellipticity;
double      MajorWidth;
double      MinorWidth;
double      MeanWidth;
double      PeakFluencePower;
int         BufferSize;
int         iShutterSetting;
double      sigCentroid[2];
double      IsoXInclusionRegionRadius_um;
double      IsoYInclusionRegionRadius_um;
double      Sigma4Ellip;
double      Sigma4EllipAngle;
double      IsoXWidth_um;
double      IsoYWidth_um;
double      ShutterSetting;
double      BaselineStd;
double      Gamma;
double      MajorWidth_dXX_WinCamD;
double      MinorWidth_dXX_WinCamD;
double      dXX_WinCamD;
double      A_dXX_WinCamD;
double      P_dXX_WinCamD;
double      IXX_WinCamD;
double      Theta_XX_WinCamD;
double      GaussianFit;
double      ImageTemp_C;
double      basic_Centroid[8];
int         Busy;
int         Minimum;
int         NumberAveraged;
int         UsedInAverage;
int         WasFullResolution;
double      PowerFactor;
char        PowerLabel[20];
double      CorrectPower;
double      InitialResult;
double      PowerInDB;
int         UseOldPowerData;
int         LogSaved;
int         MinLevel;
int         AdcPeak;
int         WasLogged;
int         Camera;
time_t      CaptureTime;
int         GammaDone;
int         Was_TwoD_Ssan;
double      PeakToAverage;
double      Ewidth_WinCamD;
int         Was_WinCamDiv;
int         SatPixels;
double      FPS;

```

```

double      EffectiveExposure;
double      PowerInCentroidTarget;
double      PlateauUniformity;
int         PixelIntensity;
double      CameraGain;
int         MatrixIndex;
double      PowerShutterSetting;
int         IsM2Data;
double      UcmM2Zlocation;
double      UcmM2SlitToLens;
double      UcmM2LensToCameraFace;
double      UcmM2LensFocalLength;
double      UcmM2WaveLength;
int         M2Data;
int         ConnectionType;
int         AdcMinimum;
double      LD;
double      ZoDelta;
double      MFactor;
int         CameraType;
int         AdcAverage;
int         PeakFound;
int         iBaseline;
int         uFIR_Gain;
int         CTE_State;
int         MeasurePeak;
int         FullResolution;
double      PowerInInclusionRegion;
int         HyperCalGood;
int         IlluminatedPixels;
int         AdcOffset;
int         Temp1;
int         ShutterState;
int         XSampleRate;
int         LineLaserCaptureWidth;
unsigned    int      AssembleTime_ms_lowByte;
unsigned    int      AssembleTime_ms_highByte;
int         IncludedGlobalWarning;
int         IntSpare;
double      TotalPower;
int         CentroidType;
double      pCentroid[2];
double      pGeoCentroid[2];
double      pPeakCentroid[2];
int         NewData;
int         ExtraLine;
BYTE        wcData[1];
}          WC_IMAGE_DATA ;

```

## Notes

1. The raw data from each frame begins at wcData [1]. Each pixel is stored as an unsigned 2-byte word and the values range from 0-65536. This does not mean the data from the camera is 16-bit. For cameras whose bitness is less than

- 16, the data is bit-shifted to fill the full range.
2. A `CPoint` and `time_t` are both 8 bytes.
3. There are some cases where Windows adds 0-padding between values. This happens for example, when there is a single `int` followed by a `double`. Windows will add 0 padding to make the `double` start 64-bits after the `int`.
4. The profiles displayed in the software are generated from the image data when the WCF file is opened. The values from the profile are not stored in the WCF file.

## Conclusion

Please contact [support@dataray.com](mailto:support@dataray.com) with any questions.