

Reading data from WCF files

The purpose of this document is to outline the structure of the .WCF file format. The recommended method of accessing data from a WCF file is to use the DataRay OCX to open the file and retrieve the calculated values as outlined in a number of tutorials on our website: <https://www.dataray.com/interfacing.html>. However, some customers need to be able to access the WCF files byte-wise.

File Structure

A WCF file contains 2 structures. First, there is a 5592 byte at the top of the file `WC_IMAGE_DATA_HEADER_2`. The first DWORD in this header contains the characters DRI. which should be used to check if the file is valid. After this header, there is a new structure for each frame stored in the .WCF file `WC_IMAGE_DATA`. This structure contains 944 bytes of header information before the actual image data. The image data is stored as a 1 dimensional array of 2-byte values. They are ordered by row. You will need to know the size of the data (rows x columns). These sizes can be read in from the Width and Height variables respectively.

```
typedef struct{
    DWORD Signature; //”DRI.”
    DWORD Type;
    DWORD Size;
    DWORD Images;
    DWORD ImagesSize;
    char Version[40];
    DRI_SETTINGS Settings;
} WC_IMAGE_DATA_HEADER_2;

typedef struct{
    int Signature;;
    int Type;
    int Index;
    int Beams;
    int Size;
    int Width;// Number of horizontal pixels
    int Height;// Number of vertical pixels
    int CameraUpdateNumber;
    double XpixelSize; //Pixel horizontal
    double YpixelSize; //Pixel vertical
    int Bits;//Normal = 16
    int Key;
    int Peak;
    int Xoffset;// x start offset (unused pixels)
    int Yoffset;// y start offset (unused pixels)
    int Xlimit;// imagers total number of x pixels
    int Ylimit;// imagers total number of y pixels
    int OreintationDone;
    CPoint pPeakCenter;
    double DefinedFluencePower;
    double pUserCentroid[2];
    double Centroid[2];
```

```
double   GeoCentroid [2];
double   Baseline;
double   UserCentroid [2];
double   GeoCenter [2];
double   PeakCentroid [2];
double   Orientation;
double   Ellipticity;
double   MajorWidth;
double   MinorWidth;
double   MeanWidth;
double   PeakFluencePower;
int      BufferSize;
int      iShutterSetting;
double   sigCentroid [2];
double   IsoXInclusionRegionRadius_um;
double   IsoYInclusionRegionRadius_um;
double   Sigma4Ellip;
double   Sigma4EllipAngle;
double   IsoXWidth_um;
double   IsoYWidth_um;
double   ShutterSetting;
double   BaselineStd;
double   Gamma;
double   MajorWidth_dXX_WinCamD;
double   MinorWidth_dXX_WinCamD;
double   dXX_WinCamD;
double   A_dXX_WinCamD;
double   P_dXX_WinCamD;
double   IXX_WinCamD;
double   Theta_XX_WinCamD;
double   GaussianFit;
double   ImageTemp_C;
double   basic_Centroid [8];
int      Busy;
int      Minimum;
int      NumberAveraged;
int      UsedInAverage;
int      WasFullResolution;
double   PowerFactor;
char     PowerLabel [20];
double   CorrectPower;
double   InitialResult;
double   PowerInDB;
int      UseOldPowerData;
int      LogSaved;
int      MinLevel;
int      AdcPeak;
int      WasLogged;
int      Camera;
time_t   CaptureTime;
int      GammaDone;
int      Was_TwoD_Ssan;
double   PeakToAverage;
double   Ewidth_WinCamD;
int      Was_WinCamDiv;
```

```

int          SatPixels ;
double      FPS;
double      EffectiveExposure ;
double      PowerInCentroidTarget ;
double      PlateauUniformity ;
int         PixelIntensity ;
double      CameraGain ;
int         MatrixIndex ;
double      PowerShutterSetting ;
int         IsM2Data ;
double      UcmM2Zlocation ;
double      UcmM2SlitToLense ;
double      UcmM2LenseToCameraFace ;
double      UcmM2LenseFocalLength ;
double      UcmM2Wavelength ;
int         M2Data ;
int         ConnectionType ;
int         AdcMinimum ;
double      LD ;
double      ZoDelta ;
double      MFactor ;
int         CameraType ;
int         AdcAverage ;
int         PeakFound ;
int         iBaseline ;
int         uFIR_Gain ;
int         CTE_State ;
int         MeasurePeak ;
int         FullResolution ;
double      PowerInInclusionRegion ;
int         HyperCalGood ;
int         IlluminatedPixels ;
int         AdcOffset ;
int         Temp1 ;
int         ShutterState ;
int         XSampleRate ;
int         LineLaserCaptureWidth ;
int         IntSpares [4] ;
double      TotalPower ;
int         CentroidType ;
double      pCentroid [2] ;
double      pGeoCentroid [2] ;
double      pPeakCentroid [2] ;
int         NewData ;
int         ExtraLine ;
BYTE        wcData [1] ;
} WC_IMAGE_DATA ;

```

Notes

1. The raw data from each frame begins at `wcData[1]`. Each pixel is stored as an unsigned 2-byte word and the values range from 0-65536. This does not mean the data from the camera is 16-bit. For cameras whose bitness is less than 16, the data is bit-shifted to fill the full range.
2. A `CPoint` and `time_t` are both 8 bytes.

3. There are some cases where Windows adds 0-padding between values. This happens for example, when there is a single int followed by a double. Windows will add 0 padding to make the double start 64-bits after the int.
4. The profiles displayed in the software are generated from the image data when the WCF file is opened. The values from the profile are not stored in the WCF file.

Conclusion

Please contact support@dataray.com with any questions.