


3x 8-32  
DEPTH 0.25 in



ACTIVE SENSOR AREA

**BladeCam-HR** 6.6 mm X 5.3 mm

**BladeCam-XHR** 6.5 mm X 4.9 mm

LED



C-MOUNT  
1.0"-32 UNC-2B

8x 4-40 UNC-2B  
DEPTH 0.2 in

USB 2.0 PORT

SENSOR PLANE

5.1 mm ±0.2  
(from sensor to case face)



DATARAY.COM  
+1-866-946-2263 +1-530-776-0843



THIRD ANGLE PROJECTION



TITLE:

BladeCam  
-HR, -XHR

	NAME	DATE
DRAWN	SM	20190503
APPROVED	RD	20190503
SIZE	LTR	

COPYRIGHT © 2018 BY DATARAY INC.

REV 3 SCALE: 1:1 SHEET 1 OF 1