


[2.400]
60.96

[0.800]
20.32


[0.000]
0




[0.000]
0

[0.925]
23.50

[2.650]
67.31



USB 2.0
Mini-B
Port




Units:
Millimeters
[Inch]

[0.000]
0

[0.870]
22.10

[2.605]
66.17



[0.000]
0

[0.925]
23.50

[1.325]
33.66

[1.725]
43.82

[2.650]
67.31

Plane Spacing	4XY models				3XYKE Models	
	100µm:	250µm:	500µm:	750µm:	50 µm:	100 µm:
First Plane	-100	-250	-500	-750	-50	-100
Second (Zero) Plane	0	0	0	0	0	0
Third Plane	100	250	500	750	50	100
Fourth Plane	400	1000	2000	3000	0	0

For a 250µm plane spacing, typical zero plane distance from front of profiler case is 3.25mm. Profiler specific distance varies by unit.

DATARAY.COM +1-866-946-2263 +1-530-776-0843			
THIRD ANGLE PROJECTION 		TITLE: <h2 style="text-align: center;">BeamMap 2</h2>	
	NAME	DATE	
DRAWN	SM	20190513	
APPROVED	LJ	20190513	
SIZE		LTR	
REV	2	SCALE: 1:1.5	SHEET 1 OF 1
COPYRIGHT © 2018 BY DATARAY INC.			