


- 4 wheels, each containing 3 filters & 1 empty space
- 4 fixed positions per wheel
- Non-parallel filters (inclined by 4°)
- C-mount threads on both ends
- Connecting adapters available
- Stray light fully eliminated
- Variable height of the optical axis
- Three mounting holes



DataRay Inc.

Advancing the Technology of Laser Beam Analysis

THIRD ANGLE PROJECTION		NAME	DATE	DWG TITLE:	
		DRAWN SM	170616	EAM-2 Reference	
		APPROVAL RD	170616		
COPYRIGHT © 2017 BY DATARAY INC.			REV: 1	SCALE:2:1	SHEET 1 OF 1